


85070E Dielectric Probe Kit Contents



	Part Number	Comes in	Description
1	85070-20037	85070E-030 and 85070-60009	Slim Form Probe (qty 3)
2	85070-60010	85070E-050	Performance Probe
3	85070-60007	85070E-030	Slim Form Probe Holder – Connects Slim Form Probe to 10mm Holder and Probe Stand Bracket .
4	85070-60011 Set	85070E	10mm Holder and Key – connects Performance Probe and Slim Form Holder to Probe Stand Bracket .
5	0905-1605	85070E-030 and 85070-60007	O-rings (qty 3) –inserted between inner and outer parts of Slim Form Probe Holder .
6	85070-20028	85070E	Probe Stand Bracket
7	85070-60003	85070E-020	High Temperature Short
8	85070-60004	85070E-030	Slim Form Short
9	8710-2036	85070E-020	High Temperature Probe
10	1250-1743	85070E	Type N to 3.5mm Adapter
11	85070-20036	85070E-020	Conductive Elastomer Short Disk – optional use adds cushion between High Temperature Probe and High Temperature Short .
12	1250-3449	85070E-030	Connector Saver Adapter – Protects Slim Form Probe when connected to 2.4mm PSC connector.
13	11901-60004	85070E	2.4mm to 3.5mm Adapter
14	85070-60008	85070E	ECal Bracket – Connects optional ECal module to Probe Stand Bracket .
15	85070-60012	85070E-050	Performance Probe Short